


SAS4305 Analytical Chemistry for Chemical Industry (AS114105/1)



## Chapter 4


### Ultraviolet-Visible (UV-Vis) Molecular Absorption Spectrometry

References:

- Skoog, D.A., Crouch, S.R., Holler, F.J., West, D.M. (2014). Fundamentals of Analytical Chemistry, 9<sup>th</sup> edition, Brooks/Cole, Chapters 24, 25 & 26.
- Skoog, D.A., Holler, F.J., Crouch, S.R. (2018). Principles of Instrumental Analysis, 7<sup>th</sup> edition, Thomson, Chapters 13 & 14.

1

Table of Content



4.1 Introduction to molecular spectroscopy

4.2 Principle of molecular absorption


4.3 Instrumentation of UV/VIS Spectrophotometer

4.4 Different design of UV/VIS Spectrophotometer

4.5 Analysis using Ultraviolet-Visible Spectrophotometer

2

4.1 Introduction to molecular spectroscopy



In Molecular spectroscopy,


- the sample solution **absorbs electromagnetic radiation** from an appropriate source,
  - and the amount **absorbed is related to the concentration** of the analyte in the solution

**Example of Molecular spectroscopy,**

- Ultraviolet and Visible Spectrophotometries (UV/VIS)
- Infrared Spectrophotometries (IR)
- Fluorescence/Phosphorescence Spectrophotometries (Flu/Phs)

3

Electromagnetic radiation (em)



The Electromagnetic Spectrum

	Gamma rays	X-rays	Ultraviolet	Visible light	Infrared	Micro waves	Radio waves
Wavelength $\lambda$ (nm)	< 0.001	1-10	190-400	400-750	750-0.1	0.1-100	100-1000
	Short wavelength			→	Long wavelength		


In Ultraviolet-visible spectrometry, the analysis based on the use of electromagnetic radiation in the wavelength region of **190 to 800 nm**

4

Electromagnetic radiation (em)

Wave properties

It is characterized with **wavelength** and **frequency**.



5

Electromagnetic radiation (em)

General terms

<b>A</b> =	<b>Amplitude</b> the length of the electrical vector at the max. in the wave
<b>P</b> =	<b>Period</b> the time required for the passage of successive max. (or min.) through a fixed point in space
<b>f</b> =	<b>Frequency</b> (unit: Hz or s <sup>-1</sup> ) no. of oscillations of the field per second
<b>λ</b> =	<b>Wavelength</b> linear distance between successive max. and min. of a wave
<b>v = f λ</b> where <b>v</b> = velocity of light in vacuum ( 3 × 10 <sup>8</sup> m/s)	
<b>v</b> =	<b>Wavenumber</b> the number of wavelengths per unit distance, cm <sup>-1</sup> <b>Wavenumber = 1 / λ</b>

6

4.2 Principle of molecular absorption

4.2.1 Molecular absorption spectrum

4.2.2 Type of energy transitions

4.2.3 Beer-Lambert Law


7

4.2.1 Molecular absorption spectrum


- more complex** than atomic absorption spectrum

→ there are **more energy states** in molecules

Atomic absorption spectrum




Molecular absorption spectrum




8

### 4.2.1 Molecular absorption spectrum

- band absorption spectrum
- e.g. from 250 to 290 nm
- molecules absorb in a **continuous range of wavelengths**



A band of  $\lambda$  (250-290nm)



Molar absorptivity,  $\epsilon$


Wavelength, nm

CN1C=NC2=C1C(=O)N(C)C2=O  
Caffeine


### 4.2.1 Molecular absorption spectrum

Molecular absorption of UV-Vis involves

- bonding electrons**,
- the **absorption band** can be used for
  - identification of different **functional groups**



$\beta$ -carotene



Molar absorptivity,  $\epsilon$


Wavelength, nm

CN1C=NC2=C1C(=O)N(C)C2=O  
Caffeine

10

### 4.2.2 Types of energy transition

- unexcited molecule **absorb specific wavelength,  $\lambda_{\text{specific}}$**
- ground state **bonding electron** transits to a higher **excited**
  - (i) **electronic**,
  - (ii) **vibrational** or
  - (iii) **rotational** states



Absorption

Energy

Excited electronic state 2


Excited electronic state 1

Vibrational energy levels

Visible

IR


Ground electronic state 1



### 4.2.2 Types of energy transition

Energy (E) being absorbed by the unexcited molecule:

$$E = E_{\text{electronic}} + E_{\text{vibrational}} + E_{\text{rotational}}$$



Absorption


Energy

Excited electronic state 2

Excited electronic state 1

Visible


Ground electronic state 1



$E_{\text{electronic}}$


### 2. Principle of molecular absorption

Energy (E) being absorbed by the unexcited molecule:

$$E = E_{\text{electronic}} + E_{\text{vibrational}} + E_{\text{rotational}}$$


### 2. Principle of molecular absorption


Energy (E) being absorbed by the unexcited molecule:

$$E = E_{\text{electronic}} + E_{\text{vibrational}} + E_{\text{rotational}}$$


### 4.2.3 The Beer-Lambert Law

**Beer's law** states that


- the **intensity** of a beam of ray of **monochromatic light (single wavelength)** **decreases exponentially** as the **concentration** of the absorbing medium **increases**.



$P_0$  - intensity of light **entering** solution  
 $P$  - intensity of light **leaving** solution  
 $b$  - length of absorbing layer (pathlength)

15

### 4.2.3 The Beer-Lambert Law



Transmittance,  $T = P/P_0 = 10^{-k'b}$   
Absorbance,  $A = -\log T$   
Absorbance,  $A = \log (P_0/P) = k'b \rightarrow A \propto b$

**100% transmittance = zero absorbance**

$P$  = intensity of light leaving solution  
 $P_0$  = intensity of light entering solution  
 $b$  = length of absorbing layer (pathlength)  
 $k'$  = constant for a particular solution

16

### 4.2.3 The Beer-Lambert Law

- the **absorbance** of the analyte is proportional to
  - (i) the **concentration** of the analyte
  - (ii) the **path length** of the sample cell

**Absorbance,  $A = \epsilon bc$**

$\epsilon$	the molar absorptivity of the solution in <b><math>\text{cm}^{-1}\text{mol}^{-1}\text{L}</math></b>
b	path length of the sample cell in <b>cm</b>
c	concentration of the solution in <b>mol/L or M</b>

17

### 4.2.3 The Beer-Lambert Law

Absorbance,  $A = \epsilon bc = \log_{10} \left( \frac{P_o}{P} \right)$

$\epsilon$	the molar absorptivity of the solution in <b><math>\text{cm}^{-1}\text{mol}^{-1}\text{L}</math></b>
b	path length of the sample cell in <b>cm</b>
c	concentration of the solution in <b>mol/L or M</b>
$P_o$	intensity of light entering solution
P	intensity of light leaving solution

18

### 4.2.3 The Beer-Lambert Law

#### Transmittance VS Absorbance

Absorbance,  $A = \log_{10} \left( \frac{P_o}{P} \right)$

Transmittance,  $T = \frac{P}{P_o}$

$A = -\log_{10} T$

**100% transmittance = zero absorbance**

★ Calculation of

- absorbance
- molar absorptivity
- path length
- transmittance(%)

19

### 4.2.3 The Beer-Lambert Law

The molar absorptivity ( $\epsilon$ )

- is a measurement of how strongly a chemical species absorbs light **at a given wavelength**.

It is an **intrinsic property** of the species.

20

4.2.3 The Beer-Lambert Law



Example

A solution with concentration  $7.25 \times 10^{-5} \text{ M}$  has a transmittance of 44.1% when measured in a 2.10 cm cell at wavelength of 525 nm. Calculate:  
(a) the absorbance of this solution  
(b) the molar absorptivity of this solution.

Solutions

(a)  $A = -\log T = -\log 0.441 = -(-0.3554) = 0.355$   
(b) Molar absorptivity  
     $= A/bc$   
     $= 0.3554 / (2.1 \text{ cm} \times 7.25 \times 10^{-5} \text{ mol L}^{-1})$   
     $= 2.33 \times 10^3 \text{ Lmol}^{-1} \text{ cm}^{-1}$

4.2.3 The Beer-Lambert Law



Classwork

A solution placed in a 1.0 cm sample cell and 70% of light is transmitted.  
If the molar absorptivity of the solution is  $2.0 \text{ cm}^{-1}\text{g}^{-1}\text{L}$ .  
What is the concentration of the solution?

22

4.2.3 The Beer-Lambert Law



Classwork (Answer)

$b = \text{_____ cm}$   
transmitted,  $T = \text{_____ \%} = 0.7$   
molar absorptivity,  $\epsilon = \text{_____ cm}^{-1}\text{g}^{-1}\text{L}$ .

$A = -\log T = \epsilon bc$


$-\log \text{_____} = (\text{_____}) (\text{_____}) (c)$   
 $\text{_____} = (\text{_____}) (c)$

$c = \text{_____} / \text{_____}$   
 $c = \text{_____ g/L}$

4.2.3 The Beer-Lambert Law



Limitation of Beer-Lambert Law



24

4.2.3 The Beer-Lambert Law



Limitation of Beer-Lambert Law

Reference ONLY!!

At **concentration > 0.01 M**, the average distances between ions or molecules of the absorbing species are diminished to the point where each particle **affects the charge distribution** and thus the extent of absorption.

25

4.3 Instrumentation of UV/VIS Spectrophotometer




- 4.3.1 main components & block diagram
- 4.3.2 radiation sources
- 4.3.3 wavelength selector
- 4.3.4 sample
- 4.3.5 detector
- 4.3.6 readout device



26

4.3.1 Main components & block diagram



5 essential components of spectrometric instruments



27

4.3.1 Main components & block diagram



★ *block diagram of UV/VIS spectrophotometers*



5 essential components of spectrometric instruments

What is a Block diagram?

Reference ONLY!!


- usually constructed to illustrate the major components of an analytical instrument.
- In a block diagram, there is no need to draw the actual structure of the instrument component

28

### 4.3.1 Main components & block diagram


5 main components in UV/VIS spectrophotometer

1. Radiation sources – stable light source
2. A wavelength selector - select interested wavelength
3. Sample cell - One or two sample cell
4. A radiation detector - converts radiant energy to measurable electrical signal.
5. A signal processing and readout device.



29

### 4.3.2 Radiation Source



**Radiation Source** → **Wavelength selector** → **Sample** → **Detector** → **Readout device**

**Properties:**


- Sufficient power
- Stable

**Ideal light source →**

- low noise
- Long-term stability


30

### 4.3.2 Radiation Source



Two kinds of radiation sources:

Continuum source	Line source
<ul style="list-style-type: none"><li>• Produce broad wavelength</li><li>• <b>NOT specific wavelength</b></li></ul>	<ul style="list-style-type: none"><li>• Emit a few discrete lines</li><li>• <b>Relatively specific</b></li></ul>




Wavelength of different radiation sources

<https://goo.gl/images/nAqWAW>

31

### 4.3.2 Radiation Source




Two kinds of radiation sources: **Reference ONLY!!**

Continuum source		Line source
In UV region	Visible & near IR region	In UV/VIS regions
For example: <ul style="list-style-type: none"><li>• Deuterium</li><li>• Hydrogen lamp</li><li>• Gas filled arc lamps,</li><li>• Light emitting diodes</li></ul>	For example: <ul style="list-style-type: none"><li>• Tungsten filament lamp</li></ul>	Example: <ul style="list-style-type: none"><li>• Mercury and sodium vapor lamps</li></ul>

32



### 4.3.3 Wavelength selector





```
graph LR; RS[Radiation Source] --> WS((Wavelength selector)); WS --> S[Sample]; S --> D[Detector]; D --> RD[Readout device];
```

- A device to isolates a desired wavelength for measurement
- No selector is capable of producing radiation of a single wavelength.
- Instead, the output of such a device is a range of contiguous wavelengths called a band.

33


### 4.3.3 Wavelength selector






<https://goo.gl/images/L2TUKi>

Interference filter



<https://goo.gl/images/f91qHg>




<https://goo.gl/images/AqgCqv>

Monochromator

34


### 4.3.3 Wavelength selector



1. **Filters**
  - a) Interference filter
  - b) Absorption filter
2. **Monochromator** → most common
  - a) Prism
  - b) Reflection grating

35

### 4.3.3 Wavelength selector



1. **Filters**
  - Provide low resolution wavelength selection (i.e. bandwidth is larger)
  - Wavelength cannot be change continuously
  - For quantitative analysis

**Example**

- a) **Interference filter:** for UV, visible & IR range
- b) **Absorption filter:** for visible range

36

### 4.3.3 Wavelength selector

#### 1. Filters

##### a) Interference filter

<https://goo.gl/images/rafDj4>

##### b) Absorption filter

<https://goo.gl/images/Uy3UDS>

37

### 4.3.3 Wavelength selector

#### 2. Monochromator

- Provide higher resolution wavelength selection (i.e. bandwidth is narrow)
- **Wavelength can be change** continuously
- For **qualitative and quantitative** analysis

Example

a) Prism

b) Reflection grating

38

### 4.3.3 Wavelength selector

#### 2. Monochromator:

Two basic designs:

##### (a) Prism

<https://goo.gl/images/ayiiyL>

##### (b) Reflection grating


<https://goo.gl/images/aBeSfq>

39

### 4.3.3 Wavelength selector

#### 2. Monochromator

a) Prism



40

### 4.3.3 Wavelength selector

**2. Monochromator**

**b) Reflection grating**

Source

Entrance slit

Concave mirrors

Reflection grating

Exit slit

Focal plane

Sample

41

### 4.3.3 Wavelength selector

**2. Monochromator:**

- A monochromator consists of lenses or mirror to focus the radiation,

Components of Monochromators:

- Entrance slit
- Collimating lens or mirror
- prisms / reflection grating
- Focal plane
- Exit slit

Source

Entrance slit

Collimating lens

Prism

Focusing lens

Focal plane

Exit slit

Sample

### 4.3.3 Wavelength selector

**1. Entrance slit:**

**Reference ONLY!!**

- protect from dust or corrosive fumes
- restrict unwanted radiation
- control the spectral bandwidth of radiation emitted from the monochromator,

Source

Entrance slit

Collimating lens

Prism

Focusing lens

Focal plane

Exit slit

Sample

### 4.3.3 Wavelength selector

**2. Collimating lens or mirror**

**Reference ONLY!!**

- produces a parallel beam of radiation

Source

Entrance slit

Collimating lens

Prism

Focusing lens

Focal plane


Exit slit


Sample

4.3.3 Wavelength selector

3. Prism or reflection grating

- disperses the radiation to different wavelengths
- by rotating the prism / reflection grating






45


4.3.3 Wavelength selector

4. Focal plane

Reference ONLY!!

- reforms the image and focuses it on a planar surface on the exit slit






46


4.3.3 Wavelength selector

5. Exit slit

Reference ONLY!!

- protect from dust or corrosive fumes
- restrict unwanted radiation
- control the spectral bandwidth of radiation emitted from the monochromator





47

4.3.3 Wavelength selector

Monochromator

a) Prism

- Prism works satisfactory in the UV/VIS region.
- By rotation of prism,
  - different wavelengths can pass the slit.

Diagram of the prism

<https://goo.gl/images/Z3he2D>

48

### 4.3.3 Wavelength selector

#### Monochromator

b) **Reflection grating**

- cheaper
- provides **better wavelength separation**
- grating consists of closely spaced grooves
- reflection on the surface due to a thin layer of **Al, Au or Pt** is coated

UV-Vis from 300 to 2000 grooves/mm  
infrared from 10 to 200 grooves/mm

Diagram of the reflecting grating

<https://goo.gl/images/k4n8bB>

49

### 4.3.4 Sample

**Sample Cell**

**Sample compartment**

reference cell

sample cell

50

### 4.3.4 Sample

#### Different Sample holder

51

### 4.3.4 Sample

#### Sample holder: Cells or cuvettes

- Materials should be transparent in the spectral region of interest
- **UV region:**
  - **Quartz or fused silica**
- **visible region:**
  - **Plastic cell**

<https://goo.gl/images/4Ay9GU>

Open-top normal with lid

52

### 4.3.4 Sample

**Blank measurement:**

- To **compensate** for the effects of absorption by the **container, reflection**, and also any **scattering effects**

The diagram illustrates the components of a sample holder. Incident radiation  $P_0$  enters from the left. Transmitted radiation  $P$  exits to the right. Various losses are indicated: reflection losses at the top and bottom interfaces, scattering losses within the solution, and reflection losses at the bottom interface. The incident beam is labeled  $P_i$  and the emergent beam is labeled  $P_e$ .

### 4.3.4 Sample

**Blank measurement:**

- To compensate the effect power beam of **analyte solution** is compared with **solvent** only

The schematic diagram shows the flow of light through a UV-Vis spectrometer. It starts with a Source (1), followed by a Wavelength selector (2), then a Sample (3) which is highlighted with a red box, then a Detector (4), and finally a Signal processor and readout (5). A red box around the Sample component contains the following text:

**Measure:**

- 1. Blank (solvent)
- 2. Standard (standard + solvent)
- 3. Sample (sample + solvent)

### 4.3.5 Detector

The flowchart shows the sequence of components in a UV-Vis spectrometer: Radiation Source → Wavelength selector → Sample → Detector → Readout device. The Detector component is highlighted with a red cloud.

- Detection of **photons**

★ Different Types of photon detectors:

3.5.1 Single channel

- (i) photo tube
- (ii) photomultiplier tube (PMT)

3.5.2 Multichannel

- (i) diode array detector(DAD)

### 4.3.5 Detector


The flowchart shows the sequence of components in a UV-Vis spectrometer: Radiation Source → Wavelength selector → Sample → Detector → Readout device. The Detector component is highlighted with a red cloud.

- Convert radiant energy into electrical signal

An ideal detector

- High sensitivity
- High signal to noise ratio
- Stabile and fast response
- Zero output signal in the absence of illumination

### 4.3.5 Detector



1. Single channel:
  - measure **single wavelength** at a time

Example:

- (i) Phototube
- (ii) photomultiplier tube (PMT)


2. Multichannel
  - Measure **different the wavelengths** simultaneously

Example:

- (i) diode array detector (DAD)

57

### 4.3.5 Detector




1. Single channel:
  - (i) Photo tube
    - An **active surface (cathode)**
    - **radiation (photon beam)**
    - causes **emission of electrons**

Diagram of the photo tube

<https://goo.gl/images/v2vVtQ>

58

### 4.3.5 Detector




1. Single channel:
  - (ii) photomultiplier tube (PMT)
    - Advantages over the ordinary phototube

Photomultiplier tube animation  
<https://www.youtube.com/watch?v=f61eMq4Wg4w>

Diagram of the photomultiplier tube (PMT)

<https://goo.gl/images/Nqo3nD>

### 4.3.5 Detector



Principle of photomultiplier (PMT)


- Photo cathode surface is similar to the surfaces of the phototubes
- Emit electrons when exposed to radiation
- Contains **additional electrode** called **dynodes**

Photomultiplier tube animation  
<https://www.youtube.com/watch?v=f61eMq4Wg4w>

Diagram of the photomultiplier tube (PMT)

<https://goo.gl/images/Nqo3nD>

### 4.3.5 Detector




#### Principle of photomultiplier (PMT)

- **Dynodes 1** is maintained at a **positive potential** higher than cathode
- Upon striking the dynodes, **each photoelectron** causes emission of **several additional electrons**
- then, **electrons** are accelerated toward **dynode 2**,
  - which is also at higher positive potential than dynode 1

Photomultiplier tube animation  
<https://www.youtube.com/watch?v=f61eMq4Wg4w>

Diagram of the photomultiplier tube (PMT)  
<https://goo.gl/images/Nqo3nD>

### 4.3.5 Detector




#### Principle of photomultiplier (PMT)

- this process has been **repeated** through dynode 1 to dynode 9
- $10^6$  to  $10^7$  electrons have been formed for each incident photon
- Cascade is finally collected at the anode
- The resulting current is **amplified** and **measured**

Photomultiplier tube animation  
<https://www.youtube.com/watch?v=f61eMq4Wg4w>

Diagram of the photomultiplier tube (PMT)  
<https://goo.gl/images/Nqo3nD>

### 4.3.5 Detector



#### 2. Multichannel


- Measure **different wavelengths simultaneously**
- Full spectrum can be generated
- Can determine more than one chemical compound simultaneously

Example:

(i) diode array detector (DAD)


63

### 4.3.6 Readout device



Radiation Source → Wavelength selector → Sample → Detector → Readout device

- electronic device that amplifies the electrical signal from the detector
- Perform mathematical operation



64



4.4 Different design of UV-VIS Spectrophotometer

4.4.1 Single-beam design

4.4.2 Double beam design

- In space (using *beam splitter*)
- In time (using *sector mirror*)

65

4.4 Different design of UV-VIS Spectrophotometer

4.4.1 Single-beam design

- with **alternate** measurement of
  - blank
  - standards
  - samples
- usually not for spectra scanning.
- Cost – cheaper than double beam design

66

4.4 Different design of UV-VIS Spectrophotometer

4.4.1 Single-beam design

1. Blank

2. Standard

3. Sample

67

4.4 Different design of UV-VIS Spectrophotometer

4.4.1 Single-beam design

Blank measurement:

- To compensate the effect power beam of **analyte solution** is compared with **solvent** only

```
graph LR; A[Radiation Source] --> B[Wavelength selector]; B --> C((Sample)); C --> D[Detector]; D --> E[Readout device]
```

Measure:

- 1. Blank (solvent)
- 2. Standard (standard + solvent)
- 3. Sample (sample + solvent)

68

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

- measurement of **sample and blank simultaneously**

69

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

★ Advantages of double beam design:

- compensate the **fluctuations** in the **radiation source**
- reduce the **fluctuations / drift** of the **detector**
- reduce the effect of **wavelength instability** of the spectrophotometer

70

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

★ in space (using beam splitter)

71

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

★ in space (using beam splitter)

- two beams are formed by a V-shaped mirror called a **beam splitter**.
- One beam** passes through the **reference solution** to a photodetector and
- the **second beam** passess simultaneously through the **sample** to a second matched photodetector.

72

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

★ in time (using sector mirror)

Readout

Reference cell

Sample cell

Grid mirror

Photo-detector

Amplifier

Source  $h\nu$

Filter or monochromator

Sector mirror

Motor

Front view

Transparent

Mirror

Sector Mirror

Transparent

73

4.4 Different design of UV-VIS Spectrophotometer

4.4.2 Double beam design

★ in time (using sector mirror)

- the **beams** are **separated in time**
- by **rotating sector mirror** that
- directs the entire beam through the reference cell
- and then through the sample cell.

Sector Mirror

Transparent

74

4.5 Analysis using UV-VIS Spectrophotometers

4.5.1 Qualitative analysis using UV/Vis

4.5.2 Quantitation analysis using UV/Vis

4.5.3 Application of UV/Vis

4.5.4 More about quantitative analysis

75

4.5.1 Qualitative analysis using UV/VIS

- changing the wavelength
- e.g. from 350 – 650nm
  - measure the absorbance
- UV spectrum (x-axis – wavelength, y-axis absorbance)**
- Compare the UV/Vis spectrum of sample/ standard/wavelength maximum

source

Wavelength selector

Change wavelength 350, 355, 360 ... 650nm

Sample

detector

Spectrum

Absorbance

Wavelength

$\lambda_{max}$

19

4.5.1 Qualitative analysis using UV/VIS



Limitation in qualitative analysis using UV-Vis

- number of absorption maxima and minima are relatively few
- difficult to identification the unknown compound

77

4.5.2 Quantitative analysis using UV/VIS



- ★ 1. Select the **wavelength maximum  $\lambda_{max}$**   
i.e. the wavelength with the maximum absorbance
2. **Fixed** the wavelength at  $\lambda_{max}$
3. **Measure the absorbance** of blank, calibration standards & sample
4. determine the concentration of sample



78


4.5.2 Quantitative analysis using UV/VIS



5. **Plot calibration curve**

- **X-axis concentration & Y-axis absorbance**

6. Determination the concentration of the sample



79


4.5.3 Application of UV/VIS



★ **Ultra-violet region**

- Mainly for analytes with **conjugated double bonds** or with **aromatic rings**
- molecules containing **unsaturated functional group**
- **Chromophore**

For example:



80

4.5.3 Application of UV/VIS



Visible region


- Mainly for analytes with colored
- Or derivatisation with coloring agent

81

4.5.4 More about quantitative analysis



- **Analysis of mixtures of absorbing substances**
- Total absorbance of a solution at any given wavelength is equal to the sum of the absorbance of the individual components in the solution



82

4.5.4 More about quantitative analysis




- **Analysis of mixtures of absorbing substances**

At  $\lambda_1$ :

$$A_{total} = A_M + A_N$$

$$A_M = \epsilon_{M_{\lambda_1}} bc_M$$

$$A_N = \epsilon_{N_{\lambda_1}} bc_N$$



83

4.5.4 More about quantitative analysis




- **Analysis of mixtures of absorbing substances**

At  $\lambda_2$ :

$$A_{total} = A_M + A_N$$

$$A_M = \epsilon_{M_{\lambda_2}} bc_M$$

$$A_N = \epsilon_{N_{\lambda_2}} bc_N$$



84

4.5.4 More about quantitative analysis

- Analysis of mixtures of absorbing substances

At  $\lambda_1$ :

$$A_{total_1} = \epsilon_{M_1}bc_M + \epsilon_{N_1}bc_N$$

At  $\lambda_2$ :

$$A_{total_2} = \epsilon_{M_2}bc_M + \epsilon_{N_2}bc_N$$

By measuring  $A_1$  and  $A_2$  and known  $\epsilon_{M_1}$ ,  $b$ ,  $\epsilon_{N_1}$ ,  $\epsilon_{M_2}$ ,  $\epsilon_{N_2}$ ,  $c_N$  and  $c_M$  can be found

85

4.5.4 More about quantitative analysis

★ **Example:**

The absorbance of solution X and solution Y is determined using UV/Vis spectrophotometer.

The results is shown in below:

Solution	Concentration (M)	Sample Cell length (cm)	Absorbance	
			400 nm	500 nm
X	$2 \times 10^{-4}$ M	1.0 cm	0.800	0.760
Y	$4.25 \times 10^{-5}$ M	1.0 cm	0.34	0.085

a) Calculate the molar absorptivity ( $\epsilon$ ) of solution X and solution Y at 400 nm and 500 nm respectively.

4.5.4 More about quantitative analysis

★ **Example (Answer):**

Using  $A = \epsilon b c \rightarrow$  Molar absorptivity,  $\epsilon = A / (b \times c)$

For **solution X** (at 400nm), molar absorptivity ( $\epsilon$ )

$$= 0.8 / (1.0 \text{ cm} \times 2 \times 10^{-4} \text{ M}) = 4000 \text{ M}^{-1}\text{cm}^{-1}$$

For **solution X** (at 500nm), molar absorptivity ( $\epsilon$ )

$$= 0.760 / (1.0 \text{ cm} \times 2 \times 10^{-4} \text{ M}) = 3800 \text{ M}^{-1}\text{cm}^{-1}$$

For **solution Y** (at 400nm), molar absorptivity ( $\epsilon$ )

$$= 0.34 / (1.0 \text{ cm} \times 4.25 \times 10^{-5} \text{ M}) = 8000 \text{ M}^{-1}\text{cm}^{-1}$$

For **solution Y** (at 500nm), molar absorptivity ( $\epsilon$ )

$$= 0.085 / (1.0 \text{ cm} \times 4.25 \times 10^{-5} \text{ M}) = 2000 \text{ M}^{-1}\text{cm}^{-1}$$

4.5.4 More about quantitative analysis

★ **Example:**

Solution X & Y then mixture together and analyze using UV/VIS spectrophotometer.

The results is shown in below:

Mixture of solution X & solution Y	Sample Cell length (cm)	Absorbance	
		400 nm	500 nm
	2.5 cm	0.680	0.520

b) Calculate the concentration (in M) of solution X and solution Y in the mixture.

	Absorbance	Sample Cell length (cm)	molar absorptivity ( $\epsilon$ )	
			Solution X	Solution Y
400nm	0.680	2.5 cm	4000	8000
500nm	0.520	2.5 cm	3800	2000

At 400nm,  $4000(2.5) C_x + 8000 (2.5) C_y = 0.680$   
 $4000 C_x + 8000 C_y = 0.680 / 2.5$   
 $4000 C_x + 8000 C_y = 0.272$

$$\begin{aligned} \text{At 500nm,} \quad & 3800(2.5) C_x + 2000(2.5) C_y = 0.520 \\ & 3800 C_x + 2000 C_y = 0.208 \end{aligned}$$

A diagram showing a beam of white light entering a triangular prism from the left. The light is dispersed into a spectrum of colors (violet, blue, green, yellow, orange, red) as it exits the prism to the right. A red dashed line indicates the path of the incident light.

$$C_v = 9 \times 10^{-6} \text{ M}$$

∴ concentration of solution X in the mixture =  $5 \times 10^{-5}$  M and concentration of solution Y in the mixture =  $9 \times 10^{-6}$  M

A diagram showing a beam of white light entering a triangular prism from the left. The light is refracted and dispersed into its constituent colors, forming a spectrum. A red circle highlights the entry point of the light, and a red arrow points from the spectrum to a small inset image showing a 3D representation of the spectrum as a rainbow-colored block.